

# WELCH ALLYN 767 SERIES

## WELCH ALLYN 767 MOBILE AND WALL ANEROIDS

The Welch Allyn 767 Series not only great looks but clinical-grade accuracy too.

- **7670-01 Wall Aneroid with Durable One-Piece Adult Cuff**
- **7670-03 Mobile Aneroid with Durable One-Piece Adult Cuff and Standard Base**



### Features:

- Certified accurate to +/- 3 mm Hg
- Mercury-free and latex-free for safety
- Ruby-jeweled movement for long life
- Laser engraved dial face for accuracy
- Large, easy to read dial face (6.25"/15 cm)
- Built-in cuff storage basket
- Mobile stand has adjustable height and 360° swivel for clear viewing
- Lifetime Calibration Warranty