

# BRONZE LINE DURASHOCK



## THE FIRST MAJOR INNOVATION IN BLOOD PRESSURE DEVICES IN 100 YEARS

The Welch Allyn Durashock range of hand-held devices provides users with a major benefit - a shock-resistant gauge which can fall 0.9m onto a hard surface and still remain accurate without requiring re-calibration! The units come complete with a standard adult cuff but will support a wide range of other sizes for maximum versatility.

### FEATURES

- Shock resistant movement
- Thumbscrew air release valve
- 5 year warranty on the gauge
- One tube friction fit

