

SILVER LINE DURASHOCK



The Welch Allyn Durashock range of hand-held devices provides users with a major benefit - a shock-resistant gauge which can fall 0.9m onto a hard surface and still remain accurate without requiring re-calibration! The units come complete with a standard adult cuff but will support a wide range of other sizes for maximum versatility.

The DS55 is available in four colour schemes - blue, yellow, red and black.

FEATURES

- Shock resistant for accuracy and lower cost over time
- Laser-engraved dial face for accuracy
- Thumbscrew air release valve
- 10 year warranty on the gauge
- Premium Aesthetics

