

The Welch Allyn Green Series™ 300 General Exam Light



1 LED



Easy Power Connector



Touchless On/Off



Easy Assembly

WelchAllyn

Advancing Frontline Care™



Green Series Medical Exam Lights

Whiter, Brighter, and Greener

Utilizing white LEDs, the GS 300 General Exam Light provides a whiter, brighter light than halogen and uses less energy consumption since LEDs never require changing.

The Welch Allyn GS 300 General Exam Light is ideal for use during general examinations, foreign body removal and minor suturing in the following environments:

- PRIMARY CARE FACILITIES
- EMERGENCY ROOMS
- OUTPATIENT FACILITIES

# of LEDs	1 White LED
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	20,000 lux (1,858 footcandles)
LED Life	50,000 hours
Working Distance	12 in (30.5 cm) to 20 in (50 cm)
Spot Size at Working Distance	5 in (13 cm) to 6 in (15 cm)
Electrical	Input: 100-240 V~ 50 Hz-60 Hz, 0.35 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Table/Wall Mount: 5.65 lbs (2.56 kg) Total Mobile Stand: 12.75 lbs (5.78 kg) Control Box: Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) Luminaire: Height: 7.48 in (19 cm) Width: 3.54 in (9 cm) Depth: 3.15 in (8 cm) Mobile Stand: Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) Total Wall Mount: 33.27 in (84.5 cm) Total Mobile Stand: 58.66 in (149.0 cm) Table/Wall and Mobile Stand to Luminaire Center: 31.1 in (79 cm)
Mounting Options	Table/Wall Mount Mobile Stand
Control Box	Touchless On/Off Control
Accessory Options	NA
Warranty	2 years

MOUNTING OPTIONS AVAILABLE:



Table/Wall Mount



Mobile Stand

To place an order, contact our Customer Service Department at 1.800.535.6663

