

The Welch Allyn Green Series™ 600 Minor Procedure Light



Disposable Sheath
Available—Use #52630



Easy Power
Connector



3 LEDs



Touchless On/Off
with Dimming Option

Easy Assembly

Green Series Medical Exam Lights



WelchAllyn®

Advancing Frontline Care™



Whiter, Brighter, and Greener

Utilizing white LEDs, the GS 600 Minor Procedure Light provides a whiter, brighter light than halogen and uses less energy consumption since LEDs never require changing.

The Welch Allyn GS 600 is ideal for use in these minor surgical procedures:

- SOFT TISSUE BIOPSY
- SKIN LESION EXCISION
- SUPERFICIAL WOUND REPAIR
- FOREIGN BODY REMOVAL
- VASECTOMY
- POST OP CARE
- BUNIONECTOMY

# of LEDs	3 White LEDs
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	50,000 lux (4,645 footcandles)
LED Life	50,000 hours
Working Distance	12 in (30.5 cm) to 20 in (50 cm)
Spot Size at Working Distance	6 in (15 cm) to 7 in (18 cm)
Electrical	Input: 100-240 V~ 50 Hz-60 Hz, 0.65 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Table/Wall Mount: 6.45 lbs (2.93 kg) Total Mobile Stand: 17.35 lbs (7.87 kg) Control Box: Height: 10.83 in (27.5 cm) Width: 3.94 in (10.0 cm) Depth: 3.15 in (8.0 cm) Luminaire: Height: 13.58 in (34.5 cm) Width: 7.56 in (19.2 cm) Depth: 3.39 in (8.6 cm) Mobile Stand: Height: 24.02 in (61 cm) Diameter: 20.87 in (53 cm) Total Wall Mount: 33.27 in (84.5 cm) Total Mobile Stand: 58.66 in (149.0 cm) Table/Wall and Mobile Stand to Luminaire Center: 31.1 in (79 cm)
Mounting Options	Table/Wall Mount Mobile Stand
Control Box	Touchless On/Off Control with Dimming
Accessory Options	Disposable Sheath #52630
Warranty	3 years



MOUNTING OPTIONS AVAILABLE:



Table/Wall Mount

Mobile Stand

To place an order, contact our Customer Service Department at 1.800.535.6663

