


The Welch Allyn Green Series™ 900 Procedure Light



6 LEDs



Touchless On/Off
with Dimming Option

Disposable Sheath
Available—Use #52630

Green Series Medical Exam Lights



WelchAllyn®

Advancing Frontline Care™



Whiter, Brighter, and Greener

Utilizing white LEDs, the GS 900 provides a whiter, brighter light than halogen and uses less energy consumption since LEDs never require changing.

The Welch Allyn GS 900 is ideal for use in the following environments:

- AMBULATORY SURGERY CENTERS
- HOSPITALS—ER AND ICU PERSONNEL
- OB/GYN
- PLASTIC SURGERY
- MINOR SURGERY OUTPATIENT FACILITIES



# of LEDs	6 White LEDs
Color Temperature	5,500° K
Light Intensity @ Typical Working Distance	80,000 lux (7,432 footcandles)
LED Life	50,000 hours
Working Distance	24 in (60.9 cm) to 40 in (100 cm)
Spot Size at Working Distance	9 in (23 cm) to 11 in (28 cm)
Electrical	Input: 100-240 V~ 50 Hz-60 Hz, 1.00 A max
Power Cord	#18 AWG, 3-wire grounded cord set, 8' (2.5 m) long
Physical Dimensions	Weight: Total Wall Mount: 17.90 lbs (8.12 kg) Total Mobile Stand: 44.85 lbs (20.34 kg) Total Ceiling Mount: 24.50 lbs (11.11 kg) Luminaire: Height: 13.74 in (34.9 cm) Width: 18.03 in (45.8 cm) Depth: 4.02 in (10.2 cm) Wall to Luminaire Center: 68.9 in (175 cm) Wall Mount Max. Up/Down: 37.4 in (95 cm) Mobile Stand: Height: 68.9 in (175 cm) Mobile Stand Base: Length: 21.97 in (55.8 cm) Width: 22.24 in (56.5 cm) Height: 6.5 in (16.5 cm) Mobile Stand Max. Up/Down: 37.4 in (95 cm) Mobile Stand to Luminaire Center: 37.4 in (95 cm) Ceiling to Luminaire Center: 68.9 in (175 cm) Ceiling Mount Max. Up/Down: 37.4 in (95 cm)
Mounting Options	Wall Mount Mobile Stand Ceiling Mount
Control Box	Touchless On/Off Control with Dimming
Accessory Options	Disposable Sheath #52630
Warranty	3 years

MOUNTING OPTIONS AVAILABLE:



To place an order, contact our Customer Service Department at 1.800.535.6663



WelchAllyn
Advancing Frontline Care™

