

Welch Allyn Maxi-Stabil 3 Sphygmomanometer



it's precise



it's robust



it's light



it's easy



WelchAllyn[®]

Advancing Frontline Care[™]

WELCH ALLYN MAXI-STABIL 3 SPHYGMOMANOMETER

Tabletop model

The Maxi-Stabil 3 sphygmomanometer is a highly versatile gauge concept, part of the Technic Series from Welch Allyn, and is designed for blood pressure measurement in a variety of professional settings.

This precision unit can be used on any tabletop or shelf. While particularly suitable for various hospital departments and clinics, it is equally at home in a physician's office or nursing home. Double-shell casing made of ABS impact polymer for shock protection.

- > Pressostabil® technology protects the internal pressure mechanism from over-inflation
- > Certified ± 3 mmHg accuracy over the full range
- > Progressive air-release valve for smooth pressure intervals
- > Easy-to-read high-visibility graduations and needle pointer
- > Velcro cuff standard, with wide range of accessory cuffs
- > Also available in latex-free design
- > Three-year warranty for movement

Maxi-Stabil 3 with Velcro cuff table-top model

- | | |
|---|-------------|
| with cuff holder | 59-34-189 |
| with cuff holder and latex-free Velcro cuff | 59-34-189LF |
| with cuff basket | 59-35-189 |
| with cuff basket and latex-free Velcro cuff | 59-35-189LF |



ALSO AVAILABLE AS:

Maxi-Stabil 3 Wall-Mounted model



- | | |
|---|-------------|
| Maxi-Stabil 3 with Velcro cuff wall-mounted model | 59-30-189 |
| Maxi-Stabil 3 with latex-free Velcro cuff | 59-30-189LF |

Maxi-Stabil 3 Rail-Mounted model



- | | |
|---|-------------|
| Maxi-Stabil 3 with Velcro cuff rail-mounted model | 59-31-189 |
| Maxi-Stabil 3 with latex-free cuff | 59-31-189LF |

Maxi-Stabil 3 Mobile Stand model



- | | |
|---|-------------|
| Maxi-Stabil 3 with Velcro cuff and stand model | 59-32-189 |
| Maxi-Stabil 3 with latex-free Velcro cuff and stand model | 59-32-189LF |

Welch Allyn GmbH & Co. KG
 Zollerstrasse 2-4
 D- 72417 Jungingen, Germany
 Tel. +49 (0) 7477-92 71-0
 Fax: +49 (0) 7477-92 71-90
 info@welchallyn.de, www.welchallyn.de