




Welch Allyn MacroView™ Otoscope

The MacroView Otoscope delivers almost twice the field of view at higher magnification compared with standard otoscopes, making it easier to visualise the auditory canal and tympanic membrane. This instrument is a significant advance in hand-held otoscopy, providing an improved view for enhanced learning and decision making.


- A nearly complete view of the tympanic membrane with approximately twice the field of view and 30% greater magnification than traditional otoscopes
- Better clarity and definition of landmarks
- The ability to adjust focus for variable ear canal length or farsighted eyes
- Tip grip for secure fastening and ejection of ear specula
- Compatible with all existing Welch Allyn 3.5 V power sources

23820 MacroView Otoscope with Throat Illuminator (Head only)

23810 MacroView Otoscope (Head only)



The MacroView Otoscope gives you a larger, clearer field of view, making it easier to diagnose conditions like otitis media.



WelchAllyn®

Advancing Frontline Care™




Welch Allyn PanOptic™ Ophthalmoscope

The revolutionary PanOptic provides easy entry into the eye, together with a wider field of view to more easily observe conditions such as hypertension, diabetic retinopathy and papilloedema.


- Fast, easy entry into undilated pupils
- 5 X larger view of the fundus vs. standard ophthalmoscopes in an undilated eye
- 25° field of view vs. the standard 5° field of view
- Greater working distance improves comfort for both practitioner and patient
- Helps detect corneal abrasions with cobalt blue filter and magnifying lens (11820)
- Compatible with all existing Welch Allyn 3.5 V power sources

11820 PanOptic Ophthalmoscope with Cobalt Blue Filter and Add-on Corneal Viewing Lens (Head only)

11810 PanOptic Ophthalmoscope (Head only)



Standard
Ophthalmoscope



PanOptic
Ophthalmoscope

5X Larger View of the Fundus



71902 3.5 V Lithium Ion Handle with
Charger and Transformer

WELCH ALLYN SOUTH AFRICA (PTY) LTD
394 Surrey Ave., Ferndale, Randburg, 2194, South Africa
P O Box 809, Fontainebleau, South Africa, 2032
Tel. +27 (0) 11 777 7555 fax +27 (0) 11 777 7556

WelchAllyn®

Advancing Frontline Care™