

The Welch Allyn Pocket Junior Otoscope

A Small Scope that Makes a Big Difference



The Welch Allyn Pocket Junior Otoscope is a valuable diagnostic tool for any office practice. Built with the same technology as Welch Allyn's other high-quality products, it provides the illumination needed to help you diagnose at an affordable price. And with its attractive styling and convenient size, it will quickly become a big part of your practice.

Benefits

- Lightweight and compact
- Halogen illumination provides brighter, whiter, longer lasting light
- Fibre optic light transmission
- Wide angle magnifying lens
- Insufflation port for pneumatic otoscopy
- Standard AA batteries with interchangeable heads for otoscope and ophthalmoscope
- Handle artwork customization available (minimum order 5000 units)

Ordering Information

22840	Pocket Junior Otoscope complete with AA handle and batteries
22841	Pocket Junior Otoscope complete with AA handle, batteries and case
03400	Replacement lamp
21504	Otoscope insufflator bulb
52432	Disposable Specula 2.5mm
52434	Disposable Specula 4mm

Presenting the Welch Allyn Pocket Junior Ophthalmoscope

Ordering Information

12850	Pocket Junior Ophthalmoscope complete with AA handle and batteries
12851	Pocket Junior Ophthalmoscope complete with AA handle, batteries and case
03900	Replacement lamp



Benefits

- Lightweight and compact
- Halogen illumination provides brighter, whiter, longer lasting light
- Sealed optics keep instrument clean and dust free
- Aperture dial offers six selections
- Glass condensing lens
- Includes red-free and cobalt blue filter
- Standard AA batteries with interchangeable heads for otoscope and ophthalmoscope
- Handle artwork customization available (minimum order 5000 units)

WelchAllyn®

Advancing Frontline Care™

The Welch Allyn Pocket Junior Set – A winning combination



Put the power, convenience, and versatility of a Welch Allyn Pocket Junior Set including an otoscope and ophthalmoscope to work for you. This set will provide the valuable diagnostic tools you need at an affordable price, and it comes in an attractive, yet durable case that's easy to carry.

95001 Pocket Junior Complete Set - Otoscope with 25 specula, Ophthalmoscope, 2 handles, batteries and case

Specifications for Welch Allyn Pocket Junior Otoscopes and Ophthalmoscopes



	Otoscope	Ophthalmoscope
illumination	Halogen	Halogen
Transmission	Fibre optic to distal tip	N/A
Batteries	2 AA	2 AA
Viewing lens	Wide angle and removable	Glass condensing lens
Weight with batteries	132g (4.65 oz.)	129.5g (4.56 oz.)
Height	16.33 cm (6.43 in)	16.02 cm (6.31 in)
Diopter Range	N/A	-25 to +22
Handle material	All metal construction	All metal construction
Insufflation port	Yes	N/A
Aperture Positions	N/A	6, including Slit and Fixation
Sealed	N/A	Yes
Filters	N/A	Cobalt Blue and Red-Free



Welch Allyn European Customer Services	Tel. (353) 46 90 67790	Fax (353) 46 90 67754
Welch Allyn UK Ltd	Tel. (44) 0207 365 6780	Fax (44) 0207 365 9694
Welch Allyn France	Tel. (33) 1 6009 3366	Fax (33) 1 6009 6797
Welch Allyn Germany	Tel. (49) 7477 92710	Fax (49) 7477 927190
Welch Allyn Italia	Tel. (39) 02 6699 291	Fax (39) 02 667 13599
Welch Allyn Iberia	Tel. (34) 91 3863047	Fax (34) 91 3861029
Welch Allyn USA	Tel. (315) 685-4560	Fax (315) 685-3361

WelchAllyn[®]

Advancing Frontline Care™