



**The 767 series features:**

- **Certified accuracy to  $\pm 3$  mm Hg.**
- **Lifetime calibration warranty.**
- **Jeweled movement for long life.**
- **Laser-engraved dial face for accuracy.**
- **Unique FlexiPort™ blood pressure cuff.**
- **Mercury-free and latex-free for safety.**
- **Large, easy-to-read dial face (15 cm).**

**Part Number & Product Description**

**767 Desktop Aneroid**

**7670-16      767 Desktop Model with Reusable One-Piece Adult Cuff**